



AXAS

(Analog X-ray Acquisition System)



Operating Manual



KETEK

Creative Detector Solutions

Custom specifications

AXAS no.00041....

Radiation entrance window:

Beryllium / AP3.3

Snout length:

50 mm / 100 mm

Detector area:

5 mm² / 10 mm²

Detector type:

SDD 5 - 145500

Serial number:

SDD 01849

Guaranteed energy resolution:

143 - 147 eV

Shaping amplifier:

yes / no

Shaping time:

150 ns / 250 ns / 500 ns / 1 µs

External power supply:

yes / no

Vacuum tightness:

yes / no

Upper energy limit:

10 keV / 20 keV / 30 keV

Other special features:

.....additional shaping time plug-in boards.
.....0.5, 0.25, 0.15 µs.
.....
.....
.....
.....

In case of questions or warranty please contact:

KETEK GmbH
Gustav-Heinemann-Ring 125
D-81739 Munich

Mail: support@ketek.net
Tel.: ++ 49 89 67346770
Fax.: ++ 49 89 67346777

Caution

- This document contains data valid only for the AXAS system you purchased it with. Please do not use it as a reference to operate any other devices you purchased earlier or later on from our company.
- Do not drop or cause mechanical shock to the AXAS, especially if you purchased a system with an open detector module.
- Avoid any kind of high forces and/or torques during mounting of the AXAS, especially if you purchased a system with an open detector module.
- Do not operate the AXAS at a temperature (measured at the steel electronics case) above 35 °C. Apply additional cooling to the steel case if necessary (contact to heat sink, gills, PC air cooler etc.). The performance of the system will noticeably worsen above this temperature.
- Do never change the settings of the sealed potentiometers inside the electronics case. Every system is adjusted carefully before delivery. We will not take any warranty for systems damaged by manipulation of these potentiometers or adjust the system a second time for free, if you did not contact us first.
- Do not touch or apply any mechanical force or shock to the radiation entrance window (Beryllium or Ap3.3) integrated at the front end of the AXAS. The integrity of each window is checked carefully before delivery and no warranty will be taken if it is damaged by the customer. Moreover breaking the window during operation of the system will most probably destroy the whole detector.
- Take care to connect the output contacts and the test input properly as it is described in this manual later on.
- Take care to follow the instructions regarding the power input to the system given in this manual. Otherwise the system may be destroyed. KETEK recommends a power supply described in this manual later on.
- Do never operate a windowless detector in a moist ambient (air is a moist ambient). It will most probably be destroyed due to the formation of water/ice on the detector surface.
- Do not touch or apply mechanical force or shock to the detector chip or mounting components if you purchased a windowless detector system. These parts are very fragile and are easily damaged.
- It is recommended to store a windowless detector under dry nitrogen atmosphere to avoid moisture and contaminations.
- Do not expose the detector to a high particle flux (such as ions or electrons). Damages done by incident particles will not be covered under warranty.
- Concerning radiation damage done by photons a maximal dose is guaranteed for each detector. In case of malfunction it is due to the customer to bring evidence that this dose has not been applied.



* The picture on the title page shows the AXAS detector system (left) and the optional power supply (right) recommended by KETEK.

Operating instructions

The AXAS-detector system should be kind of plug and play. That means you should not have to change anything on the system itself. If this though is necessary contact us first in any case! After you unpacked the different components please connect them to each other and the components you supply by yourself regarding to the connection diagram inside this manual.

Packing list:

You should have received the following components:

- AXAS detector system, including detector, preamplifier, shaping amplifier (removed on customers demand), temperature control and internal power supply
- External ± 12 V power supply if purchased.
- Additional shaping time boards if purchased.
- This manual. (Invoice will be sent with separate mail!)

Dimensions and SUB-D connector:

A technical drawing inside this manual gives the dimensions of the AXAS system, the pin definition of the integrated male SUB-D9 connector, and the definitions of the other input and output connectors in detail. Please never connect any pin, input or output in a different way or you may destroy the detector. Please note that the drawing corresponds to an AXAS system with a 100 mm snout. So this dimension may vary if you purchased a system with a different snout length.

Shaped output / radiation entrance window:

The BNC / Shaper OUT output provides an amplified, shaped signal (see example inside this manual / typical rise and fall time curves are attached too) within a voltage range from 0 to the value you specified (maximal 10 V, input impedance 1 mega ohm). It can be used as input signal to a standard ADC. The detector is calibrated in a way that the mentioned voltage spread covers a spread of x-ray energies from 0 to the upper energy limit you specified. The theoretical lower limit of the detector is about 200 eV but due to the entrance window you purchased the practical lower limit may be higher (e.g. for Beryllium about 1 keV). For more detailed information on the transmission of your entrance window please have a look at the attached transmission curves.

Pre-amplified output / test pulse input:


If you did not purchase our shaper the BNC / Shaper OUT does not provide any signal. In any case a pre-amplified signal is provided by BNC / Preamp. OUT. BNC / Test Input allows the user to connect a pulse generator to the system and test the output by applying different test pulses.

Pile-up rejection:

The pins 3 and 5 of the SUB-D9 connector each provide a TTL signal as input for an ADC with integrated pile-up rejection (not_PUR is inverted with respect to PUR / see attached oscilloscope screen shot).

AXAS 00041

SSD01849



AXAS Datasheet



KETEK GmbH
Halbleiter- und
Reinraumtechnik

Every detector system produced by KETEK GmbH is tested electrically and spectroscopically. Testing parameters and results achieved in the final qualification are reported below.
A leakage test (MIL-STD-883, Method 1014.10) was performed and passed.

AXAS-Serial-NO: 00041
 AXAS-Batch-NO: 2002.02
 TO8-Serial-NO: 01849
 Class: SDD5-145500

Spectroscopic results

FWHM: 145 eV @ 1 μ s @ -20 °C
 (chip temperature)

P/B ratio: 382

Count rate [cps]: 996

TO8 housing

Cap Type:	MOXTEK DuraBeryllium 8 um
Collimator:	Dual collimator
Collimator support:	Covar + Gold Layer

Entrance window: 2.4 mm

Auxiliary Coll. dia.:

1. Sealing
 The TO8-housing is closed hermetically. A leakage test MIL-STD-883 was performed. However every mechanical stress on the detector module can damage the module entrance window or the PIN's introducing a leakage in the housing.


2. Radiation hardness

Radiation hardness of the silicon drift detector has been investigated using synchrotron radiation.
 With a total dose of 1012 absorbed photons on the SDD active area, typical variation of the energy resolution in terms of FWHM are less than 2% and peak/background ratio variation are less than 10%.
 Radiation hardness is energy dependent. The analyses made represents the worst case, assuming incoming x-rays in an energy range between 10 and 15 keV, where the damage on the integrated electronics is maximized. Outside this energy range is the detector radiation harder.

The spectroscopic results can be achieved, when the temperature of the AXAS housing is +20°C.

Our technical staff is at your disposal for any kind of questions
 mail to: support@keteck.net
<http://www.keteck.net>

Technical drawing of an AXAS system with a 100 mm snout



Ambient temperature:

Please note that the detector performance depends on the ambient temperature you provide. The metal housing of the AXAS system acts as a heat sink for the cooling mechanism of the detector. The housing should not have a temperature above 35 °C during operation. Anyway it is optimal if the housing stays at room temperature or somewhat below. So if possible attach an additional heat sink or apply active cooling to the housing.

Detector temperature:

You can monitor the current temperature of the SDD detector integrated into your AXAS system by reading the voltage given at the two temperature monitor output connectors (see technical drawing). This voltage is calibrated to deliver 3.00 V at 300 K and will in-/decrease 0.01 V per K (e.g. 273 K = 2.73 V or 312 K = 3.12 V).

Resolution (at high count rates) / shaping time:

Each detector has been tested carefully at the KETEK laboratory at Munich/Germany to redeem the guaranteed specifications. The guaranteed energy resolution was checked with reference to a measurement using the Mn K_α line of a Fe55 x-ray source and a shaping time of 1 μs. Please note that the best resolution you will observe depends on the shaping time and the count rate you use. For more detail have a look at the attached diagram giving the resolution of an example AXAS system in dependence of count rate and shaping time.

Peak shift at high count rates:

For high count rates a peak shift in the recorded spectra may occur depending on the applied count rate. The attached “AXAS – Peak Shift” diagram gives an example of this peak shift. Please note that the example we attached represents a kind of worst case scenario. The actual peak shift is different for every detector and should not exceed the given example by much.

Windowless detector:


If you purchased a windowless detector never operate it during venting or at any moist atmosphere (air is a moist atmosphere) or it will most probably be destroyed.

Check the connection diagram:

If your system does not operate properly, please check all the connections referring to the connection diagram. If everything is connected the right way, please contact us before you open the AXAS and change any internal settings. **YOU MAY DESTROY THE DETECTOR!!!**




Connection diagram




*It is recommended to use Keteks shaping amplifier integrated in the AXAS system


Example spectrum of a Fe55 radiation source $\sigma(\text{Mn K}_\alpha) = 141 \text{ eV FWHM}$




AXAS - FWHM



AXAS - Peak Shift




Double pulse and PUR TTL pulse




This screen shot shows a typical double (pile up) pulse (upper curve) and the subsequent TTL pulse produced by the AXAS system (lower curve) as input to an optional pile up rejection feature of your ADC. Please note, that the original output is inverted with respect to this picture.

Amplified, shaped output of the AXAS system




This oscilloscope screen shot shows a typical shaped and amplified output signal of the AXAS system (BNC / Shaper OUT). A Fe⁵⁵ x-ray source was used for irradiation of the detector. The Mn K_α and K_β lines are clearly separated. Please note that the amplitude of the signal depends on the adjusted gain voltage range and the x-ray energy.

Rise time of the pre-amplified signal




This picture shows the typical rise time (139.1 ns in this case) of a pre-amplified signal (BNC / Preampl. OUT). Again a Fe⁵⁵ x-ray source was used to irradiate the detector. Please note that the amplitude and the total rise time of the signal depend on the x-ray energy and that the original output is inverted with respect to this picture.

Fall time of the pre-amplified signal




This picture shows the typical fall time (98.44 μ s in this case) of a pre-amplified signal (BNC / Preamp. OUT). Again a Fe⁵⁵ x-ray source was used to irradiate the detector. Please note that the amplitude and the total fall time of the signal depend on the x-ray energy and that the original output is inverted with respect to this picture.

Transmission of a AP3.3 window



Transmission of a duraberyllium window




1.008	H	Atomic Mass	
Li	6.941 4 K α 0.054	Be K α 0.109	9.012 K α 0.4461 K β 0.395
Na	22.99 12 K α 1.041	Mg K α 1.254	24.31
K	39.10 20 K α 3.314 K β 3.590	Ca K α 3.692 K β 4.013 L α 0.341	40.08 21 K α 4.091 K β 4.461 L α 0.395
Rb	85.47 38 K α 13.40 K β 14.96 L α 1.694	Sr K α 14.17 K β 15.84 L α 1.807	87.62 39 K α 14.96 K β 16.74 L α 1.923
Cs	132.9 56 K α 30.97 K β 34.97 L α 4.287 L β 4.620	Ba K α 32.19 K β 36.38 L α 4.466 L β 4.828	137.3 57 K α 33.44 K β 37.80 L α 4.651 L β 5.042
Fr	(223) 88 L α 12.03 L β 14.77 M α 2.743	Ra L α 12.34 L β 15.24 M α 2.825	(226) 89 Ac L α 12.65 L β 15.71 M α 2.910



KETEK

Creative Detector Solutions





2 He		4.003																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																															
5	B K α 0.183	6	C K α 0.277	7	N K α 0.392	8	O K α 0.525	9	F K α 0.677	10	Ne K α 0.849	11	20.18																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																				
13	Al K α 1.487	14	Si K α 1.740 K β 1.836	15	P K α 2.014 K β 2.139	16	S K α 2.308 K β 2.464	17	Cl K α 2.622 K β 2.816	18	Ar K α 2.958 K β 3.191	19	39.95																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																				
9	Sc K α 6.41 24 K β 6.490 L α 0.637	10	Ti K α 4.511 K β 4.932 L α 0.452	11	V K α 4.952 K β 5.427 L α 0.511	12	Cr K α 5.115 K β 5.947 L α 0.573	13	Mn K α 5.899 K β 6.490 L α 0.637	14	Fe K α 6.404 K β 7.058 L α 0.705	15	Co K α 7.478 K β 8.265 L α 0.776	16	Ni K α 8.639 K β 9.905 L α 0.930	17	Cu K α 8.048 K β 9.056 L α 0.852	18	Zn K α 8.639 K β 9.577 L α 1.012	19	Ga K α 9.252 K β 10.26 L α 1.098	20	Ge K α 10.54 K β 11.73 L α 1.282	21	As K α 11.22 K β 12.50 L α 1.379	22	Se K α 11.92 K β 13.29 L α 1.480	23	Br K α 12.65 K β 14.11 L α 1.586	24	Kr K α 12.65 K β 14.11 L α 1.586	25	83.80																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																
7	Ru K α 19.38 43 K β 21.66 L α 2.697 L β 2.834	8	Pd K α 20.22 K β 22.22 L α 2.839 L β 2.990	9	Rh K α 19.38 K β 21.66 L α 2.559 L β 2.633	10	Pt K α 21.18 K β 23.82 L α 2.839 L β 3.151	11	Ag K α 22.16 K β 24.94 L α 3.134 L β 3.317	12	Os K α 23.17 K β 25.28 L α 3.444 L β 3.663	13	In K α 24.21 K β 27.28 L α 3.769 L β 4.030	14	Sn K α 26.36 K β 28.74 L α 3.938 L β 4.221	15	Tl K α 27.47 K β 31.00 L α 4.110 L β 4.423	16	Te K α 28.61 K β 32.29 L α 4.221	17	I K α 29.78 K β 33.62 L α 4.410 L β 4.423	18	Xe K α 30.78 K β 34.92 L α 4.513 L β 4.622	19		20		21		22		23		24		25		26		27		28		29		30		31		32		33		34		35		36		37		38		39		40		41		42		43		44		45		46		47		48		49		50		51		52		53		54		55		56		57		58		59		60		61		62		63		64		65		66		67		68		69		70		71		72		73		74		75		76		77		78		79		80		81		82		83		84		85		86		87		88		89		90		91		92		93		94		95		96		97		98		99		100		101		102		103		104		105		106		107		108		109		110		111		112		113		114		115		116		117		118		119		120		121		122		123		124		125		126		127		128		129		130		131		132		133		134		135		136		137		138		139		140		141		142		143		144		145		146		147		148		149		150		151		152		153		154		155		156		157		158		159		160		161		162		163		164		165		166		167		168		169		170		171		172		173		174		175		176		177		178		179		180		181		182		183		184		185		186		187		188		189		190		191		192		193		194		195		196		197		198		199		200		201		202		203		204		205		206		207		208		209		210		211		212		213		214		215		216		217		218		219		220		221		222		223		224		225		226		227		228		229		230		231		232		233		234		235		236		237		238		239		240		241		242		243		244		245		246		247		248		249		250		251		252		253		254		255		256		257		258		259		260		261		262		263		264		265		266		267		268		269		270		271		272		273		274		275		276		277		278		279		280		281		282		283		284		285		286		287		288		289		290		291		292		293		294		295		296		297		298		299		300		301		302		303		304		305		306		307		308		309		310		311		312		313		314		315		316		317		318		319		320		321		322		323		324		325		326		327		328		329		330		331		332		333		334		335		336		337		338		339		340		341		342		343		344		345		346		347		348		349		350		351		352		353		354		355		356		357		358		359		360		361		362		363		364		365		366		367		368		369		370		371		372		373		374		375		376		377		378		379		380		381		382		383		384		385		386		387		388		389		390		391		392		393		394		395		396		397		398		399		400		401		402		403		404		405		406		407		408		409		410		411		412		413		414		415		416		417		418		419		420		421		422		423		424		425		426		427		428		429		430		431		432		433		434		435		436		437		438		439		440		441		442		443		444		445		446		447		448		449		450		451		452		453		454		455		456		457		458		459		460		461		462		463		464		465		466		467		468		469		470		471		472		473		474		475		476		477		478		479		480		481		482		483		484		485		486		487		488		489		490		491		492		493		494		495		496		497		498		499		500		501		502		503		504		505		506		507		508		509		510		511		512		513		514		515		516		517		518		519		520		521		522		523		524		525		526		527		528		529		530		531		532		533		534		535		536		537		538		539		540		541		542		543		544		545		546		547		548		549		550		551		552		553		554		555		556		557		558		559		560		561		562		563		564		565		566		567		568		569		570		571		572		573		574		575		576		577		578		579		580		581		582		583		584		585		586		587		588		589		590		591		592		593		594		595		596		597		598		599		600		601		602		603		604		605		606		607	</

AXAS

Setup and Operation


Each unit of drift detector and AXAS board is shipped with supply voltages adjusted for optimum energy resolution. **It is not recommended to change any of the voltage settings.**

In case the correct voltage settings have been lost or a detector shall be used together with a differently adjusted electronics board adjust the potentiometers as showed in the picture below. Every short circuit which can arise during this operation may destroy the silicon drift detector.



KETEK

Backside Connectors



Creative Detector Solutions


page 2 of 4

Changing the shaping time constant on AXAS system


CAUTION

The AXAS PCBs contain static-sensitive components. Handle the boards at a static-free work area.

- disconnect electrical connections
- open the case
- remove shaper unit as shown in figure



- change left and right PCB
- the shaping time constant is written on the solder side of the PCBs



- reinsert the shaping module



Detector Voltage and Other Trimming

The 6 trimmers behind the test points (the PCB detector connector) adjust the same first 6 voltages that are listed in the test point table below. Rotating the trimmer screws clockwise increases all the voltages. The Back voltage is derived from the Ring 18 voltage so the R18 voltage has to be trimmed first, provided that the Back voltage is initially decreased. After trimming the R18 voltage, the back voltage is increased to the desired value.

test points:

1	2	3	4	5	6	7	8	9	10	11	12	13	14
R18	Back	Reset	TDGR	R1	IGR	TD-n	TD-p	+12V	GND	GND	P+	-9V	FFD

As comparatively high ohm resistors are used in the detector voltages filter networks, it is a good idea to use a voltmeter of higher than 10Mohm input impedance. Standard DMM can be used if a 1:100 voltage divider is realized by connecting of a 1Gohm 1% resistor serially to the active DMM terminal. This is most convenient to do at the DMM side. Doing so, the voltages can be measured without errors due to the DMM input impedance.

Rotating the temperature control set trimmer clockwise cools the detector. The system is delivered with maximum cooling. The relation between the voltage of the temperature monitor output on the back side and the temperature is 10mVper degree K. For instance 300K (27degrees C) are read as 3.000V. An offset of this temperature monitor can be adjusted with the "temp. monitor calibration" trimmer.

Fine gain control trimmer clockwise rotation increases the gain of the shaping amplifier. However, it should be kept in mind that the Coarse Gain (fixed) is about 0.25 V/keV (0.5, 0.33 or 0.17 V/keV). With a maximum output voltage of 5 V the shaper gives signal for up to 20 keV (10, 15 or 30 keV).

The trimmer on the preamplifier board can be used for fine P/Z adjustment.

