

mar research

The **mar345** Software Guide

By Dr. Claudio Klein



X-ray Research G.m.b.H.
Segeberger Chaussee 34
D-22850 Norderstedt
Germany

Tel.: (+49) (40) 529-884-0
Fax: (+49) (40) 529-884-20
Email: info@marresearch.com

The *mar345* Software Guide

Version 1.2
April 4, 2001

Written by Dr. Claudio Klein

Copyright 1997–2001 X–Ray Research G.m.b.H.

All Rights Reserved.

This document contains proprietary and confidential information of X–ray Research G.m.b.H.. The contents of this document may not be disclosed to third parties, copied or duplicated in any form, in whole or in part, without the prior written consent of X–ray Research G.m.b.H..

Possession, use, or copying of the software described in this publication is authorized only pursuant to a valid written license from X–ray Research G.m.b.H.

marresearch and the *marresearch* logo are registered trademarks of X–ray Research G.m.b.H..

Contents

1	Distribution Media	1
2	Compilation Notes	1
3	Environment	1
4	Directory Structure	2
5	Description of Programs	3
6	Documentation	4
7	Installation on SGI	5
8	Installation on DEC Unix	7
9	Installation on Linux	9
10	Running mar345	11

1. Distribution Media

Media	OS	Suggested directory	Mount CD-Rom with
CD-Rom	IRIX 5.x,6.x	/usr/people/mar345	mount -t iso9660 -r /dev/rdisk/dksXdYvol /CDROM
	DEC Unix 5.0	/usr/users/mar345	mount -t iso9660 -r -o rrip /dev/rzYc /mnt
	Linux, 2.0.x	/home/mar345	mount -t iso9660 -r /dev/cdrom /mnt/cdrom

X denotes the SCSI-controller (e.g.0) and Y the SCSI-id (e.g. 4).

Most executable programs are also available as gzip compressed files via anonymous FTP from ftp.EMBL-Hamburg.de, directory pub/marx/mar345/.... Due to their size (180 MB), scanner calibration files are not available per FTP. Please contact X-ray Research for assistance.

2. Compilation Notes

OS	Compilers used	Motif version used	
IRIX	C-7.1	Motif 1.2.4	X11R6
DEC Unix	C-5.2	Motif 1.2.4	X11R5
Linux	GCC-2.7.2	Motif 2.1.x	X11R6

3. Environment

The following logicals assignments must be set to run certain programs:

Variable	Description	Used by
MARHOME	Master directory of software distribution.	–
MARLOGDIR	Directory where log files will be created.	<i>mar345, mar345xf, scan345, marHKL</i>
MARMANDIR	Directory with man pages.	–
MARHELPDIR	Directory with online help files.	<i>mar345, marView</i>
MARTABLEDIR	Directory with scanner specific tables.	<i>marTools</i>
MAR_SCANNER_NO	Scanner serial number.	<i>mar345, mar345xf, scan345</i>

4. Directory Structure

The **mar345** software distribution directory (\$MARHOME) contains the following subdirectories:

Subdirectory	Contents
bin/sgi bin/osf bin/linux	Binary executables for corresponding platforms.
man/1 man/cat1 man/man1 man/ps man/doc man/help	Unformatted man pages for selected programs. Formatted man pages (SGI). Compressed unformatted man pages (Linux & DEC Unix). Postscript files of formatted man pages, ready for printing. ASCII text of formatted man pages, ready for online read (more). Online help files for GUI's (<i>mar345</i> , <i>marView</i> , <i>marTools</i>).
log log/log log/lp log/spy log/sets	Log-files for programs <i>mar345</i> and <i>scan345</i> . Up to 99 versions of mar.log files. Up to 99 versions of mar.lp files (statistical output). Up to 99 versions of mar.spy files (native controller messages). Data collection strategy template files for <i>mar345</i> .
src tables	Source code of selected programs. Scanner specific calibration and configuration files.
Optional: marflm marhkl marxds	Latest <i>marFLM</i> distribution. Latest <i>marHKL</i> distribution. Latest <i>marXDS</i> distribution.

5. Description of Programs

Name	man page	Priority	Description
<u>Graphical User Interfaces:</u>			
mar345	mar345	A	GUI for data collection and display.
marstart	–	A	Works together with program <i>mar345</i> and must be available.
marView	marView	A	GUI for data display.
marTools	marTools	A	GUI for image file manipulation: format conversion, etc.
tkmarcv	tkmarcv	A	GUI for image format conversion. Relies on Tcl/Tk package (not installed in a standard SGI installation).
marHKL	–	B	GUI for HKL processing package (DENZO, xdisp).
marFLM	–	B	GUI for MOSFLM processing package.
marXDS	–	B	GUI for XDS processing package (XDS, xscale).
<u>Hardware related programs:</u>			
mar345xf	mar345xf	C	Standalone transformation program for spiral images.
marsim	marsim	D	Scanner simulator. Works together with programs <i>mar345</i> or <i>scan345</i> . This is for performance tests and/or debugging.
scan345	scan345	D	Non-GUI data collection program: no display, no goniometer functions, just scans and spiral to Cartesian conversion.
modnb	–	D	Modifies header of calibration files.
swapnb	–	D	Swaps bytes in calibration file. Useful when replacing the data collection computer by one with a different byte-order.
<u>Format conversion programs:</u>			
marcv	marcv	A	Conversion of <i>mar345</i> formats into other formats. Also used by tkmarcv (see above).
spiral(un)pack	spiralpack	C	(De-)Compression of spiral images as produced by <i>mar345</i> .
<u>Others:</u>			
catmar	catmar	B	Dumps image header (<i>mar345</i> and <i>mar300</i> formats, also for calibration files).
marstats	–	C	Dumps average intensity and sigmas of images.
<u>Data processing related programs for use within <i>marHKL</i>, <i>marFLM</i> and <i>marXDS</i>:</u>			
marPeaks	marPeaks	B	Spot search program.
marIndex	marIndex	B	Autoindexing program (based on marPeaks).
marPredict	marPredict	B	Pattern prediction program (based on marIndex).
marStrategy	marStrategy	B	Program to calculate an optimal data collection strategy.

Priority codes:

- A Essential
- B Helpful, installation recommended
- C Not essential, may be removed
- D Needed only in special situations.

6. Documentation

The documentation can be found in directory \$MARHOME/man. Several formats are available:

Directory	Description
man/1	Unformatted man pages for selected programs.
man/cat1	Formatted man pages (SGI).
man/man1	Compressed unformatted man pages (Linux & DEC Unix).
man/ps	Postscript files of formatted man pages, ready for printing.
man/doc	ASCII text of formatted man pages, ready for online read (<i>more</i>).
man/help	Online help files for GUI's (<i>mar345</i> , <i>marView</i> , <i>marTools</i>).

The following man pages are available:

Name	Description
<u>Programs:</u>	
mar345	Documentation for program <i>mar345</i> .
marView	Documentation for program <i>marView</i> .
marTools	Documentation for program <i>marTools</i> .
tkmarcv	Documentation for program <i>tkmarcv</i> .
mar345xf	Documentation for program <i>mar345xf</i> .
scan345	Documentation for program <i>scan345</i> .
marsim	Documentation for program <i>marsim</i> .
marcv	Documentation for program <i>marcv</i> .
spiralpack	Documentation for program <i>spiralpack</i> .
catmar	Documentation for program <i>catmar</i> .
marPeaks	Documentation for program <i>marPeaks</i> .
marIndex	Documentation for program <i>marIndex</i> .
marPredict	Documentation for program <i>marPredict</i> .
marStrategy	Documentation for program <i>marStrategy</i> .
<u>Others:</u>	
mar345_formats	Documentation for <i>mar345</i> image formats.
mar300_formats	Documentation for <i>mar300</i> image formats.
mar345_config_file	Documentation for the config file entries for program <i>mar345</i> .

To view the man pages using man, the directory \$MARHOME/man must be in the man page search path. Consult the "man" man page for further details, since this varies from computer to computer.

The GUI's provide "Help"-buttons for additional online information.

When run with the "-h" command line option, usage information is provided for most of the mar programs, e.g. type:

```
marcv -h
```

7. Installation on SGI IRIX

7.1 Installing the Software (IRIX 6.5)

7.1.1 Login as user "root" and create a new user account called *mar345*.

Do the following:

- a) From the toolchest (usually in the upper left corner of the screen) choose **"System Manager"** from the **"System"** menu.
- b) In the main window of the **"System Manager"**, choose: **"Security and Access Control..."**
- c) Select **"Add a User Account"** and go through steps 1 to 10.
 - Choose **"mar345"** as **"Login name"** (suggested).
 - Choose **"Password"**, **"User-ID"** and **"Primary Group"** as desired.
 - Choose **"/usr/people/mar345"** as **"Home Directory"** (suggested).
 - Choose **"tcsh"** as **"Shell Program"** (strongly recommended).
 - As a last step, press **"OK"** to accept the settings.

7.1.2 Login as user: *mar345*

7.1.3 Installing the software from CD-ROM:

Insert the CD-ROM in the CD-ROM reader. After a couple of seconds, the automounter should mount the CD automatically as /CDROM. If not, the super-user should try:

```
mount -t iso9660 -r /dev/rdisk/dks0d4vol /CDROM
```

where 4 is the SCSI-unit number of the CD-ROM on the system, which may vary. If you are not certain, which SCSI-unit your CD-ROM drive uses, try command:

```
hinv
```

to get a "hardware inventory" of your system which will say something like:

```
CDROM: unit 4 on SCSI controller 0:
```

When successful, the contents of the CD-ROM should be copied into the login directory of the account *mar345*. To do so, as user "*mar345*" type:

```
/CDROM/mar_install
```

The installation script chooses reasonable defaults that may be accepted or modified.

It is important that the contents of the CD-Rom are really copied to the login directory of the new user since the distribution contains customized startup files (.cshrc, etc.) which should reside in the login directory.

7.1.2 Setting Up the Ethernet Connection

The scanner control program *mar345* communicates with the scanner through a Ethernet interface. To use *mar345*, the network must be configured to meet the requirements of the mar controller. The mar scanner will respond only to requests made to IP-address 192.0.2.1. The mar scanner also requires the host computer Ethernet interface to be set to address 192.0.2.2. These addresses belong to a pool of addresses that are not assigned to official networks so there should not be any conflict with the outside world.

Do the following:

7.1.2.1 Edit file /etc/hosts, e.g. type:

```
jot /etc/hosts
```

and add the following line to the end of the file:

```
192.0.2.1    mar345
192.0.2.2    gate-$HOSTNAME
```

where \$HOSTNAME is the original name of your computer, e.g. if the computer is called "mars", then write:

```
192.0.2.2    gate-mars
```

7.1.2.2 If the second network interface is not yet configured, edit file /etc/config/netif.options:

```
jot /etc/config/netif.options
```

and modify lines:

```
: if2name=
: if2addr=gate-$HOSTNAME
```

to:

```
if2name=ec1
if2addr=gate-$HOSTNAME
```

7.1.2.3 Reboot.

7.1.2.4 Check settings

The second network interface (ec1) should now be configured and you should be able to access other hosts (e.g mar345) on network 192.0.2. To check, type:

```
/usr/etc/ifconfig ec1
```

This command should come back with something like:

```
ec2: flags=c63<UP,BROADCAST,NOTRAILERS,RUNNING,FILTMULTI,MULTICAST>
    inet 192.0.2.2 netmask 0xfffff00 broadcast 192.0.2.255
```

The correct routing table can be checked using command:

```
netstat -rn
```

It should say something like:

Destination	Gateway	Netmask	Flags	Refs	Use	Interface
default	192.168.0.1		UGS	1	0	ec0
127.0.0.1	127.0.0.1		UH	5	592	lo0
192.0.2	link#1	0xfffff00	UC	0	0	ec1
192.168	link#2	0xfffff00	UC	0	0	ec0
196.168.0.3	127.0.0.1		UGHS	2	2	lo0
224	link#2	0xfffff00	UCS	0	0	ec0
255.255.255.255	192.168.0.255		UGHS	0	0	ec0

To manually configure ec1, type:

```
/usr/etc/ifconfig ec1 inet 192.0.2.2
route add net 192.0.2.0 gate-$HOSTNAME 0
```

Note: You may also use the graphical system configuraton tools ("System" -> "System Manager" -> "Network and Connectivity" -> "Network Interface Manager" -> "Ethernet ec1") to configure the second Ethernet card. In practice, the GUI's have failed so many times, that it seems much easier to do the steps 7.1.2.1 and 7.1.2.2 manually.

8. Installation on DEC Unix

8.1 Installing the Software (DEC Unix 4.0)

8.1.1 Create a new user account called *mar345*.

Type:

- a) **su** = become superuser
- b) **/usr/bin/X11/dxaccounts** = call the "Account Manager" graphical interface
- c) In the main window of the "**Account Manager**", choose: "**Add...**"
- d) In the window "**Create/Modify Local User**" fill in at least the "**User name**" and "**Shell**" press "**OK**".
As default shell, the `/bin/tcsh` is recommended. However, the `"tcsh"` is not installed by default on DEC Unix 4.0. If `"tcsh"` is unavailable (look also in `/usr/bin` and `/usr/local/bin`) choose: `/bin/csh`.
As "**Home Directory**" use of `/usr/users/mar345` is suggested. At this time you may create a password for this user account but you don't have to. All other fields do not need to be modified.
Press "**OK**" to accept the settings.
- f) Press "**Close**" to leave the "**Account Manager**".
- g) Log out.

8.1.2 Login as user: *mar345*

8.1.3 Installing the software from CD-ROM:

Insert the CD-ROM in the CD-ROM reader. Since normally there is no automounter, the mounting probably needs to be done as super-user. Type:

```
mount -t cdfs -o rrip -r /dev/rz0c /mnt
```

When successful, the contents of the CD-ROM should be copied into the login directory of the account *mar345*. To do so, as user "*mar345*" type:

```
/mnt/mar_install
```

The installation script chooses reasonable defaults that may be accepted or modified. It is important that the contents of the CD-Rom are really copied to the login directory of the new user since the distribution contains customized startup files (`.cshrc`, etc.) which should reside in the login directory.

8.2 Setting Up the Ethernet Connection

The scanner control program *mar345* communicates with the scanner through a Ethernet interface. To use *mar345*, the network must be configured to meet the requirements of the mar controller. The mar scanner will respond only to requests made to IP-address 192.0.2.1. The mar scanner also requires the host computer Ethernet interface to be set to address 192.0.2.2. These addresses belong to a pool of addresses that are not assigned to official networks so there should not be any conflict with the outside world.

Do the following:

8.2.1 Configure the second Ethernet card (tu1)

Type:

- 1.) **su** = become superuser
- 2.) **/usr/sbin/netconfig** = call the "Network Configuration" GUI
- 3.) Choose **"Ethernet Interface – tu1"**
- 4.) **Select:** **"tu1 – Interface Configuration:"** **Enable**
- 5.) Enter **"Host name:"** **marcontrol**
- 6.) Enter **"IP address:"** **192.0.2.2**
- 7.) Enter **"Network Mask:"** **255.255.255.0**
- 8.) Press: **Commit**

8.2.2 Edit file /etc/hosts, e.g. type:

vi /etc/hosts

and add the following line to the end of the file:

192.0.2.1 mar345

If you can't find an entry for name marcontrol, also add:

192.0.2.2 marcontrol

8.2.3 Close the "Network Configuration" GUI.

The second network interface (tu1) should be configured now and you should be able to access other hosts (e.g mar345) on network 192.0.2. To check, type:

ifconfig tu1

This command should come back with something like:

```
tu1: flags=c63<UP,BROADCAST,NOTRAILERS,RUNNING,MULTICAST,SIMPLEX>
    inet 192.0.2.2 netmask ffffff00 broadcast 192.0.2.255 ipmtu 1500
```

The correct routing table can be checked using command:

netstat -r

It should say something like:

Routing tables

Destination	Gateway	Flags	Refs	Use	Interface
Netmasks:					
Inet	255.255.255.0				
localhost	localhost	UH	5	592	lo0
192.0.2	marcontrol	U	1	0	tu1
193.141.161	cal	U	47	779	tu0

8.2.4 Connect the scanner to the second Ethernet card and power it up.

9. Installation on Linux (RedHat Distribution)

9.1 Installing the Software

9.1.1 Create a new user account called *mar345*.

Type:

- a) **su** = become superuser
- b) **adduser mar345** = adds user account "mar345" using defaults.
- c) Edit **/etc/passwd** and change the default shell of user "**mar345**" to **/bin/tcsh**.
Possibly, the password of account "**mar345**" has to be set.

9.1.2 Login as user: *mar345*

9.1.3 Installing the software from CD-ROM:

Insert the CD-ROM in the CD-ROM reader. If there is an automounter, the CD-Rom is going to be mounted automatically (on RedHat usually as **/mnt/cdrom**). Otherwise, on many Linux system, users are allowed to mount CD-Roms by just typing:

	mount	/mnt/cdrom	(RedHat)
or	mount	/cdrom	(SuSE)

If this doesn't work, the super-user has to do something like:

```
mount -t iso9660 -r /dev/cdrom /mnt/cdrom
```

When successful, the contents of the CD-ROM should be copied into the login directory of the account *mar345*. To do so, as user "*mar345*" type:

```
/mnt/cdrom/mar_install
```

The installation script chooses reasonable defaults that may be accepted or modified. It is important that the contents of the CD-Rom are really copied to the login directory of the new user since the distribution contains customized startup files (.cshrc, etc.) which should reside in the login directory.

9.2 Setting Up the Ethernet Connection

The scanner control program *mar345* communicates with the scanner through a Ethernet interface. To use *mar345*, the network must be configured to meet the requirements of the mar controller. The mar scanner will respond only to requests made to IP-address 192.0.2.1. The mar scanner also requires the host computer Ethernet interface to be set to address 192.0.2.2. These addresses belong to a pool of addresses that are not assigned to official networks so there should not be any conflict with the outside world.

Do the following:

9.2.1 Configure the second Ethernet card (eth1)

- Type:
- 1.) **su** = become superuser
 - 2.) **/usr/bin/netcfg** = call the "Network Configuration" GUI
 - 3.) Choose section **"Interfaces"**
 - 4.) Choose button **"Add"**
 - 5.) Choose **"Interface Type:"** **Ethernet**
 - 6.) In window **"Edit Ethernet/Bus Interfaces"**
 - Enter: **"IP":** **192.0.2.2**
 - Enter: **"Netmask":** **255.255.255.0**
 - Select: **"Activate interface at boot time"**
 - 8.) Press: **Done**
 - 9.) Confirm: **Save configuration**
 - 10.) In main window, select interface **"eth1"** and click button **"Activate"**
 - 11.) Click button **"Save"** and **"Quit"**.

For the kernel to actually initialize 2 Ethernet cards, please consult the documentation of the Linux boot load **"lilo"**, i.e. read the lilo man page or the lilo documentation in **/usr/doc/lilo***.

Note: On SuSE systems, the integration of the second Ethernet card is done via the SuSE system administration tool (yast or yast2). Program netcfg is not available.

9.2.2 Edit file **/etc/hosts**, e.g. type:

vi /etc/hosts

and add the following line to the end of the file:

192.0.2.1 mar345

If you can't find an entry for name marcontrol, also add:

192.0.2.2 marcontrol

The second network interface (eth1) should be configured now and you should be able to access other hosts (e.g mar345) on network 192.0.2. To check, type:

ifconfig eth1

This command should come back with something like:

```
eth1      Link encap:10Mbps Ethernet  HWaddr 00:80:C6:FF:EF:08
          inet addr:192.0.2.2  Bcast:192.0.2.255  Mask:255.255.255.0
          UP BROADCAST RUNNING MULTICAST  MTU:1500  Metric:1
          RX packets:0 errors:0 dropped:0 overruns:0
          TX packets:0 errors:0 dropped:0 overruns:0
          Interrupt:12 Base address:0x320
```

The correct routing table can be checked using command:

netstat -r

It should say something like:

Kernel IP routing table

Destination	Gateway	Genmask	Flags	MSS	Window	irrt Iface
193.141.161.0	*	255.255.255.0	U	1500	0	0 eth0
192.0.2.0	*	255.255.255.0	U	1500	0	0 eth1
127.0.0.0	*	255.0.0.0	U	3584	0	

9.2.3 Connect the scanner to the second Ethernet card and power it up.

10. Running mar345

10.1 Turn the scanner on. Wait until you hear the plate locking sound (click!).

10.2 Login as user mar345.

10.3 Try to talk to the scanner, i.e. type:

```
ping mar345 or /usr/etc/ping mar345 (SGI)
```

If the scanner is accessible, ping comes back with:

```
PING mar345 (192.0.2.1): 56 data bytes
64 bytes from 192.0.2.1: icmp_seq=0 ttl=255 time=1 ms
...
```

If ping comes back with:

```
ping: mar345: Unknown host
```

then, mar345 has not been inserted into file /etc/hosts (see 2.1).

If ping hangs with:

```
PING 192.0.2.1 (192.0.2.1): 56 data bytes
```

then the reason might be

- a) the interface has not been configured correctly.
- b) the scanner is not turned or not yet ready to listen.
- c) there is a problem with the Ethernet cable.
- d) there is a problem with the T-connectors.
- e) there is a problem with the scanner itself.

10.4 In a terminal window, type "go" to start the mar software (mar345).