



Argonne
NATIONAL
LABORATORY

... for a brighter future



U.S. Department
of Energy

UChicago ►
Argonne_{LLC}



**Office of
Science**

U.S. DEPARTMENT OF ENERGY

A U.S. Department of Energy laboratory
managed by UChicago Argonne, LLC

GM/CA-CAT Canted Undulator Beamlines for Macromolecular Crystallography: Current and Future Capabilities

Bob Fischetti
Associate Director, GM/CA-CAT
Biosciences Division
Argonne National Laboratory

January 24, 2007

APS Cross Cut Review of Biological Sciences

Funding and Community Served

GM/CA CAT is funded by the US National Institutes of Health's partnering institutes the National Institute of General Medical Sciences (**GM**), and the National Cancer Institute (**CA**) to build and operate a facility for macromolecular crystallography.

Beamtime is allocated in the following categories

- GM special grantees – Protein Structure Initiative (20.5%)
- CA special grantees – Structure based drug design (12.5%)
- General user facility (including DOE 25%) (42.0%)
- CAT Director, SAB members, R&D and staff research (25.0%)

GM/CA CAT Staff

Administration



Janet Smith
Director



Robert Fischetti
Associate Director



Sheila Trznadel
Admin Specialist

Crystallographic Support



Ward Smith



Nukri Sanishvili



Michael Becker



Naga Venugopalan

Computing Support



Sergey Stepanov



Oleg Makarov



Mark Hilgart

Engineering & Technical Support



Shenglan Xu



Derek Yoder



Rich Benn



Steve Corcoran

Scientific & Technical Vision

Wide energy tunability for MAD

3.5 - 35 keV (wavelength 3.5 - 0.35 Å)

Variable beam size to match to sample size

Easily varied in the range 10 - 200 µm

Small beam capability

Down to 10 µm routinely

High stability of beam position and wavelength

Small beam stable throughout a user's experiment

Wavelength tunable without beam position changing

User-friendly beamline operation

Easily understood & operated computer interface

Automation to speed sample screening

For both difficult and high-throughput projects, and for beamline setup

Not just high throughput


Larger macromolecules & complexes, small crystals, weakly diffracting crystals such as membrane proteins

GM/CA-CAT Capabilities


- Two independent, rapidly tunable, canted-undulator beamlines
 - 23-ID-B: 3.3 cm undulator, energy range: 3.5 – 20 keV
 - 23-ID-D: 3.0 cm undulator (optimized for Se and Br), energy range: 5.0 – 35 keV with monochromator crystal change
- One rapidly tunable bending magnet beamline, energy range: 3.5 – 35 keV with monochromator crystal change (being commissioned)
- Small, stable beams on ID-lines - bimorph mirrors arranged in a K-B geometry
 - Small focal size 20 x 70 microns
 - Homogeneous beam profile at focus and at sample
- ALS-style sample automounters on all 3 beamlines
- Excellent sample viewing through an on-axis microscope
- Convenient user interface: SSRL's Blu-Ice GUI adapted to EPICS control environment
- Mini-beam apparatus for small (< 10 micron) crystals
- CCD detectors: MAR 300 on ID-line, MAR 225 on BM-line

1st APS Sector Based on Hard X-ray Dual Canted Undulators

Typical Single Undulator




Dual Canted Undulator




- Independent steering: 3 RF-BPMs
Stable: 1 X-ray BPM at 22 m
Increased flexibility:
- 3.3 cm – standard UA
 - 3.0 cm – optimized for Se and Br



Theoretical Tuning Curves for 3.0-cm and 3.3-cm Undulators
Low Emittance Mode (3.0 nm-rad, 1% coupling, 7 GeV, 100 mA)



Schematic of Canted ID Beamlines



Milestones

ACCEL construction contract signed

White light delivered to FOE

White beam to ID-C

First monochromatic light ID_{in}

First diffraction data on ID_{in}

First user on ID_{in}

KB “bimorph” mirrors installed ID_{in}

First monochromatic light from dual canted undulators

First diffraction experiments on ID_{out}

First monochromatic light on BM (under going commissioning)

Sample automounter made available to users

First mini-beam data set collected on users sample

May 20, 2003

August 11, 2003

April 20, 2004

June 14, 2004

July 8, 2004

July 12, 2004

August 16, 2004



February 16, 2005

December 21, 2005


February 23, 2006

December 1, 2006

December 16, 2006




ACCEL - Double Crystal Monochromator



2nd Crystal

1st Crystal






- Constant exit height
- Long 2nd crystal to minimize tune error
- 2nd crystal translates in Bragg perpendicular direction
- 1st and 2nd crystals are both indirectly, cryogenically cooled
- Compton scatter shield around 1st crystal and 2nd crystal mount for improved thermal stability
- In-vacuum Huber 430 rotation axis – avoids rotary feedthrough
- Si (111) or Si(311) for higher energy and/or resolution

Thin web under cut

Reduces absorbed power and amount of Compton scattering produced

Monochromator Performance

- Energy resolution in good agreement with Si (111)
- High heat load tolerance
 - No broadening of rocking curve with 600 W on 1st crystal at 100 mA
 - Can accept 90% of central cone at 200 mA (FEA calculations)
- Minimal rocking curve dispersion or “detuning” of 1st and 2nd crystals vs. energy
- Calibration procedures for non-Bragg axes to stabilize beam position to less than $\pm 5 \mu\text{m}$ over the range 4 – 20 keV
- Spectroscopy-grade edge scans
- BPM after DCM for intensity feedback and YAG crystal for visualizing the beam



K-B “Bimorph” Mirrors Performance – achromatic focusing

Early results – August 2004

Focused at the sample

21 x 62 microns (FWHM)


Focused 1 m downstream of sample

26 x 91 microns (FWHM)

Energy Independence

Beam profile at 12 keV (red)



Beam profile at 18 keV (blue)



Uniform electrode voltages – set curvature



Differential electrode voltages – correct slope error

Bimorph longitudinal cross section



On-Axis Crystal Viewing Optics – YAG to visualize the beam

- Motor controlled 16X zoom
- Beam coincident with optical axis (through high numerical aperture objective lens)
- Field of view - 2.4 mm X 3.2 mm down to 0.15 mm X 0.20 mm
- Front and backlight




150 μm

200 μm

Easy alignment of sample to beam

- Center rotation axis on cross hair
- Move beam to cross hair with slide bar control of mirror pitch angle (piezo actuators)
- Mount sample and center on cross hair


Precision Omega Axis – small sphere-of-confusion



- AeroTech ABR1000 air bearing:
<1 μm wobble
- XY sample positioner:
2 – 3 micron sphere of confusion


Sample Environment

Beamstop "Flipper" Assembly



Beam stops automatically moves out of the way when hutch door opens, and repositions when the door closes.




Motorized scatter guard



Scatter guard moves out of optical path to allow visualization of the crystal. Data collection software repositions before taking data.

IO monitor after fast shutter

- Al scatter foil and PIN Diode
- Positioned after fast shutter
- Provides beam intensity measurement during exposure



Captured Top-Up


Automated sample handling

Modified ALS robot design, collaborated with LRDesign to incorporate improvements suggest by Thomas Earnest, Carl Cork and others, and to customize for our goniometry.




- Installed on all three beamlines
- Available to users on ID-lines
- Compatible with either ALS/SSRL (96 pins), or Rigaku (72 pins) pucks
- Screening takes < 3 min/crystal with “point & click” centering
- Implementing automated sample centering (hiring Sadhir Babu Pothineni from Victor Lamzin’s group (EMBL/Hamburg) - XREC
- Configuration allows manual mounting

Blu-Ice user interface



SSRL's Blu-Ice GUI adapted to EPICS control environment

Blu-Ice user interface



SSRL's Blu-Ice GUI adapted to EPICS control environment

Blu-Ice user interface

The screenshot displays the Blu-Ice GUI window titled "BLU-ICE Beamline gmca23id". The interface includes a menu bar (File, Component, Shutter, Network, View, Options, Help) and a tabbed view with "Screening" selected. Below the tabs, there is a "Cassette" dropdown set to "left: undefined" and an "Update" button. The main area contains a table with columns: Mount, Dismount, Port, ID, Protein, Comment, and Directory. The table lists 18 components, each with a "Mount" and "Dismount" button. Below the table is a log window showing several messages from January 16, 2006, at 12:43:27 and 12:43:36, including warnings and notes about client activity. At the bottom, a status bar shows "Idle", "Undulator: 8.2690 keV", "Energy: keV", "Network: Active", "Shutter: Closed", and the time "02:48:21 PM".


	Mount	Dismount	Port	ID	Protein	Comment	Directory
>	Mount	Dismount	A2	c_A2	0	0	0
1	Mount	Dismount	A3	c_A3	0	0	0
2	Mount	Dismount	A4	c_A4	0	0	0
3	Mount	Dismount	A5	c_A5	0	0	0
4	Mount	Dismount	A6	c_A6	0	0	0
5	Mount	Dismount	A7	c_A7	0	0	0
6	Mount	Dismount	A8	c_A8	0	0	0
7	Mount	Dismount	B1	c_B1	0	0	0
8	Mount	Dismount	B2	c_B2	0	0	0
9	Mount	Dismount	B3	c_B3	0	0	0
10	Mount	Dismount	B4	c_B4	0	0	0
11	Mount	Dismount	B5	c_B5	0	0	0
12	Mount	Dismount	B6	c_B6	0	0	0
13	Mount	Dismount	B7	c_B7	0	0	0
14	Mount	Dismount	B8	c_B8	0	0	0
15	Mount	Dismount	C1	c_C1	0	0	0
16	Mount	Dismount	C2	c_C2	0	0	0
17	Mount	Dismount	C3	c_C3	0	0	0
18	Mount	Dismount	C4	c_C4	0	0	0

16 Jan 2006 12:43:27 WARNING: This client is passively viewing the beamline.
16 Jan 2006 12:43:27 NOTE: This client may now issue commands.
16 Jan 2006 12:43:36 WARNING: This client is passively viewing the beamline.
16 Jan 2006 12:43:37 NOTE: This client may now issue commands.
16 Jan 2006 12:43:43 WARNING: This client is passively viewing the beamline.
16 Jan 2006 12:43:44 NOTE: This client may now issue commands.

Idle Undulator: 8.2690 keV Energy: keV Network: Active Shutter: Closed 02:48:21 PM 4 AM PM

SSRL's Blu-Ice GUI adapted to EPICS control environment



Blu-Ice user interface




SSRL's Blu-Ice GUI adapted to EPICS control environment

Ramp Up of the User Program

Total Number for all Runs



	Crystal Screened	Sets Collected	PDBs
ID-D	10979	2382	62
ID-B	2355	469	8
Total	13334	2852	70







GM/CA-CAT Accomplishments

- 26 Articles in refereed journals
 - *1-Cell*
 - *1-Nature*
 - *7-PNAS*
 - *2-Nature Chem Biol,*
 - *1-Nature Struc Mol Biol*
- 1 Cover (Structure)
- 70 Structures deposited (PDB)
- 44 Invited talks
- 2 Awards
- 3 Conference proceedings
- 3 Press releases on the GM/CA facility (C&E News, R&D Magazine, NIH News)

Large Unit Cells

Diffraction pattern from HK97 virus capsid.
Unit cell dimensions: 1010 x 1010 x 732 Å




MAR 225
S-D distance 680 mm

L. Gan, *et al.* & J. E. Johnson *Structure* 14, 1655-65 (2006)

Large Unit Cells, continued

Model fitting into electron density map of HK97 capsid (Jack Johnson/Lu Gan)

“...Figure below shows a section of electron density from the 4.2Å map of Head II. The quality is superb and bulky side chains can be visualized even at this resolution....”



Mini-beam Capability

Requirements

- ~5 μm beam
- improved beam positional stability
- accurate beam position in visual field
- effective lighting for ~5 μm colorless samples
- sample position stable to < 5 μm
- goniometer center of rotation stable to < 5 μm
- easy alignment of sample to beam

All of the above attributes must be stable through sample changes.

Advantages of mini-beam


- Raster long crystals
- Select best part of crystal
- Reduced background
- Projects that only produce small crystals

Disadvantages of mini-beam

- Diffraction is weak from small crystals
- Radiation damage is still a factor


Smaller and smaller beams: the practical challenge

Most common approach to a small beam for crystallography: overfill an aperture



- Tolerates beam positional instabilities
- Misses lots of flux
- Increased background from slits
- Not suitable for micro-focus

Our approach to a small beam for crystallography: focus to desired size




- Can provide high resolution data
- Captures lots of flux
- Good signal-to-noise
- Suitable for micro-focus
- Sensitive to beam positional instabilities


Mini-beam approach combines both: focus to small beam at the sample position and then aperture to the desired beam size (similar to ESRF ID13 and micro-diff)

- Beam defining slits are typically too far away ($> 200\ \text{mm}$)
- Mostly reduce intensity but not size of beam at sample


Mini-beam apparatus: small beam-defining pinholes



- Kinematics mount
- rapidly interchangeability
 - 15 minute alignment




5, 10, 30 μm apertures




7 μm beam (FWHM)
Viewed on YAG crystal at sample position
Beam focused at sample position
10 μm aperture

Knife Edge Scans of “Mini-beam” – 10 micron pin hole

FWHM = 6.3 microns



FWHM = 7.8 microns



Flux through pin hole: $\sim 1 \times 10^{11}$ photons/sec/100 mA

Still need to attenuate!

5 micron pin hole – 3.8 x 6.8 microns FWHM (V x H)

Visualization of small crystals

7 x 7 x 15 μm^3 lysozyme crystal in 30 micron MiTeGen loop






Image recorded with mini-beam apparatus in position

Diffraction patterns from small crystals

7x7x15 μm^3 lysozyme crystal




10x4x<4 μm^3 thaumatin crystal




Diffraction from small crystals is weak, but measurable

Refined data from micro-crystals

Region of lysozyme structure refined at 2.25 Å resolution against the data collected from a single small crystal.



Region of thaumatin structure refined at 1.94 Å resolution against the data collected from 6 small crystals and then merged.





Electron densities are plotted with 1 sigma cut-off.

Improved Signal/Noise Ratio from Same Crystal

7- μm beam
2.4 Å

75- μm beam
2.6 Å



Same
crystal &
setting

Same
reflection

Future


Improved efficiency

- Already have many automated tasks
- Developing automated alignment of optical components
- Developing automated focusing of mirrors
- Developing automated alignment of goniometry and centering of sample

Usage


- Lots more screening – small crystals in particular
- Remote access
- Opens the door for project stuck on the shelf because of small or weakly scattering crystals

Structural Biology Beamlines at the Advanced Photon Source



Percent of time available to general users

GM/CA Sector Layout



K-B “Bimorph” Mirrors Performance - off focus beam profiles



Focus at 72.3 m : 0 mm

Sample at 72.0 m : 300 mm

Guard slits at 71.8 m : 500 mm


Horizontal beam shape

Vertical beam shape




Only half of the electrodes are independent
More power supplies installed January 2007

Schematic of BM Beamline




Energy Range: 3.5 – 25 keV
 Flux: $0.5 - 1.2 \times 10^{12}$ photons/sec/100mA
 Focal size: 110 x 150 microns
 Status: commissioning



Beamline Automation

- Parallel homing: brings beamline to life in less than 10 minutes
- FPE audit records actual exposure time and dose per frame
- Beam intensity feedback stabilizes $I(0)$
- Beam positional feedback locks position
- Fast-scan scripting library written for automated alignment & fast XAS scans
- Auto-alignment of optics



Gallery of PDB deposits from IDin

70 deposits for which all or part of data were collected at 23IDD or 23IDB lines

